



DETERMINAZIONE DI PFAAs IN MATRICI AMBIENTALI: ESPERIENZE E CRITICITA'.



Acque del Chiampo s.p.a.
Servizio Idrico Integrato


CRONISTORIA ACCREDITAMENTO PFAAs in AdC: Metodiche e matrici

- 2015: EPA 537 «*Acque destinate al consumo umano*»: 14 isomeri lineari + 4 isomeri ramificati
- 2016: EPA 537 «*Acque destinate al consumo umano*»: 5 isomeri ramificati di PFOS e PFOA
- 2017: ASTM D 7979 «*Acque di scarico trattate*»: isomeri lineari
- 2018: ASTM D 7979 «*Acque superficiali*»: isomeri lineari
- 2019: EPA 537.1 «*Acque destinate al consumo umano e di falda*»: isomeri lineari + ramificati
- 2020: ASTM D 7979 «*Acque destinate al consumo umano*»: isomeri lineari



ANALISI PFAAs EFFETTUATE IN AdC A PARTIRE DA AGOSTO 2013:

- 2013: 373
- 2014: 1702
- 2015: 1704
- 2016: 2717
- 2017: 3725
- 2018: 4343
- 2019: 2974
- 2020: 1909



Acque del Chiampo s.p.a.
Servizio Idrico Integrato

METODI ADOTTATI DA AdC:

- EPA 537-2009 «Determination of selected perfluorinated alkyl acids in drinking water by solid phase extraction and liquid chromatography/tandem mass spectrometry (LC/MS/MS)»
- EPA 537.1-2018 «Determination of selected per and polyfluorinated alkyl acids in drinking water by solid phase extraction and liquid chromatography/tandem mass spectrometry (LC/MS/MS)»
- EPA 8327-2019 «Per and polyfluoroalkyl substances (PFAS) using external standard calibration and multiple reaction monitoring (MRM) liquid chromatography/tandem mass spectrometry (LC/MS/MS)»
- ASTM D 7968-17a «Determination of Polyfluorinated compounds in soil by liquid chromatography tandem mass spectrometry (LC/MS/MS)»



Acque del Chiampo s.p.a.
Servizio Idrico Integrato

METODI ADOTTATI DA AdC:

- ASTM D 7979-19 «Determination of per and polyfluoroalkyl substances in water, sludge, influent, effluent, and wastewater by liquid chromatography tandem mass spectrometry (LC/MS/MS)»
- ISO 25101-2009 «Determination of perfluorooctanesulfonate (PFOS) and perfluorooctanoate (PFOA) – Method for unfiltered samples using solid phase extraction and liquid chromatography/mass spectrometry»
- UNI CEN/TS 15968-2010 «Determinazione del perfluoro ottano sulfonato (PFOS) estraibile in articoli solidi, rivestiti e impregnati, nei liquidi e nelle schiume antincendio- Metodo per il campionamento, l'estrazione e l'analisi per mezzo di LC-QMS o LC-MS»
- RAPPORTI ISTISAN 19/7 «Metodi analitici per il controllo delle acque da destinare e destinate al consumo umano ai sensi del DL.vo 31/2001 e s.m.i. – ISS.CBA.051/052.REV00»



Acque del Chiampo s.p.a.
Servizio Idrico Integrato

MATRICI ANALIZZATE E RISPETTIVI LoQ:

- Acque destinate al consumo umano: < 5 ng/l
- Acque di scarico trattate + acque superficiali: PFOS < 5ng/l, altri< 10ng/l
- Acque di scarico non trattate: PFOS < 5ng/l, altri < 10ng/l
- Acque sotterranee/falda: < 5 ng/l
- Rifiuti liquidi/fanghi: PFOS < 20ng/l, altri < 50ng/l
- Rifiuti solidi/fanghi: < 0,5 µg/Kg
- Percolati discarica: PFOS < 5ng/l, altri < 10ng/l
- Prodotti chimici: 50 ng/l ÷ 0,5 µg/Kg
- Suoli/terreni: < 0,5 µg/Kg
- Eluato/test di cessione: PFOS < 5ng/l, altri < 10ng/l
- Pelli/tessuti: < 0,1 µg/Kg



CIRCUITI INTERLABORATORIO:

- 2013: n.d.
- 2014: n.d.
- 2015: 4 (acque potabili) di cui 3 z-score ok, 1 warning + JRC
- 2016: n.d.
- 2017: 2 (acque naturali), 1 z-score ok, 1 warning
- 2018: 2 (acque scarico/superficiali), z-score ok
- 2019: 11 (acque scarico/superficiali), z-score ok
- 2020: 6 (acque scarico/superficiali) + 8 (lisciviato), z-score tutti ok



ALTRI PFAAs IMPLEMENTATI:

- PFPA, PFHpS, PFNS, PFDS, PFDoS, PFTrA, PFTeA, PFPeS, P35DMHxS, P4MHpA, P5MHpA, P6MHpA, P55DMHxA, P1MHpS, P3MHpS, P4MHpS, P5MHpS, P6MHpS
- Gen x ≡ HFPO-DA
- C6O4
- 4:2 FTS, 6:2 FTS, 8:2 FTS (Fluorinated Telomer Sulfonic acid)



Acque del Chiampo s.p.a.
Servizio Idrico Integrato

INTERVENTI SOSTENUTI DA ADC:

Analisi acque potabili/pozzi per cittadini comuni soci

Interventi su reti e impianti

Predisposizione punti di distribuzione acqua non contaminata (fontanelle pubbliche)

Interventi di interconnessione

Completamento sistemi idrici

Estensione rete idrica per utenti non serviti da acquedotto




Investimenti attuali e futuri

INSTALLAZIONE SISTEMI FILTRAZIONE A CARBONI ATTIVI «GAC»

- 2013 (Agosto): Centro Idrico Madonna dei Prati-Brendola
- 2014 (Marzo): Impianto filtrazione carboni attivi Lonigo
- 2019: Centro Idrico Canove (predisposto per futura filtrazione GAC)
- 2020: Centrale di filtrazione di Montorso



ACIDI PERFLUOROALCHILICI (PFAAs): RENDIMENTO FILTRI A CARBONE ATTIVO



Acque del Chiampo s.p.a.
Servizio Idrico Integrato

PERFLUOROALKYL ACIDS




WHAT DID WE LEARN?



Acque del Chiampo s.p.a.
Servizio Idrico Integrato



PFOA

PERFLUOROOCTANOIC ACID



PFOS

PERFLUOROOCTANSOLFONIC ACID



Acque del Chiampo s.p.a.
Servizio Idrico Integrato



...some chemical and physical data

Abbreviation	$\log K_{AW}$	$\log K_{OW, \text{dry}}$	$\log S_W (\text{mol L}^{-1})$
Perfluorinated carboxylate acids (PFCAs)			
PFBA	-3.23	2.82	0.42
PFPeA	-2.90	3.43	-0.37
PFHxA	-2.58	4.06	-1.16
PFHpA	-2.25	4.67	-1.94
PFOA	-1.93 (-3.0)	5.30	-2.73
PFNA	-1.58	5.92	-3.55
PFDA	-1.27	6.50	-4.31
PFUnDA	-0.92	7.15	-5.13
PFDoDA	-0.58	7.77	-5.94
PFTrDA	-0.38	8.25	-6.59
PFTeDA	0.03	8.90	-7.42
Perfluorinated sulfonic acids (PFSAs)			
PFBS	-2.59	3.90	-1.00
PFHxS	-2.38	5.17	-2.24
PFOS	-1.65	6.43	-3.92




Acque del Chiampo s.p.a.
Servizio Idrico Integrato

...chromatographic retention time...







Acque del Chiampo s.p.a.
Servizio Idrico Integrato

...chromatographic retention time...



Acque del Chiampo s.p.a.
Servizio Idrico Integrato

Solid Phase Extraction (SPE):



Acque del Chiampo s.p.a.
Servizio Idrico Integrato



T_R – RECOVERY - PRECISION

Analyte	Retention time	RECOVERY	RSD%	INTERVAL CONCENTRATION
		n = 4 - 7		
PFBA	5,7	95 – 101%	4 - 10%	200 - 50000
PFBS	11	85 – 122 %	3 - 7%	200 - 50000
PFPeA	10,2	97 – 115%	3 - 8%	200 - 50000
PFHxA	15	99 – 128 %	2 - 7%	200 - 50000
PFHpA	18,8	97 – 127 %	4 - 14%	200 - 50000
PFHxS	19	92 – 129 %	2 - 13%	200 - 50000
PFOA	21,7	91 – 102 %	2 - 14%	200 - 50000
PFNA	24,1	84 – 106 %	7 - 19%	200 - 50000
PFOS	24,2	75 – 123 %	6 - 20%	200 - 50000
PFDA	26,1	76 – 125 %	10 - 20%	200 - 2000
PFUnA	27,8	76 – 129 %	3 - 14%	200 - 2000
PFDoA	29,3	97 – 121 %	5 - 13%	200 - 2000



Acque del Chiampo s.p.a.
Servizio Idrico Integrato

SURROGATES MIX

$^{13}\text{C}_4\text{-PFBA}$
 $^{13}\text{C}_4\text{-PFOA}$
 $^{13}\text{C}_4\text{-PFOS}$
 $^{13}\text{C}_3\text{-PFHxS}$
 $^{13}\text{C}_2\text{-PFDA}$

INTERNAL STANDARD MIX



$^{13}\text{C}_2\text{-PFHxA}$
 $^{13}\text{C}_3\text{-PFBA}$
 $^{13}\text{C}_8\text{-PFOA}$
 $^{13}\text{C}_8\text{-PFOS}$
 $^{13}\text{C}_3\text{-PFBS}$
 $^{13}\text{C}_2\text{-PFDoA}$



Acque del Chiampo s.p.a.
Servizio Idrico Integrato

L-PFOS: 499 > 80 (SO₃)



Voltaggi: cono = 60V; cc = 39V



**Acque del Chiampo s.p.a.
Servizio Idrico Integrato**

L-PFOA: 413 > 369

Voltaggi: cono = 16V; cc = 10V



Acque del Chiampo s.p.a.
Servizio Idrico Integrato

PERFLUOROALKYL ACIDS

**WHAT
CRITICAL ISSUES
?**




Acque del Chiampo s.p.a.
Servizio Idrico Integrato

TRANSITIONS AND VOLTAGES

MS file			
Analyte	Transition	Cone voltage (V)	Collision cell voltage (V)
PFBA	213 > 169	18	10
PFPeA	263 > 219	15	9
PFOS	499 > 80	60	39
	499 > 99	60	39
PFBS	299 > 80	54	30
	299 > 99	54	28
PFDA	513 > 219	20	16
	513 > 269	20	16
PFNA	163 > 219	20	18
	163 > 169	20	16
PFHpA	363 > 119	16	20
	363 > 169	16	16
PFHxA	313 > 119	16	17
	313 > 269	16	10
PFUnA	563 > 219	20	20
	563 > 269	20	18
PFOA	413 > 369	16	16
	413 > 169	16	10
PFHxS	399 > 80	40	40
	399 > 99	40	40
PFDoA	613 > 169	25	15
	613 > 319	25	15
	613 > 569	25	15




PFOS ISOMERS: CHROMATOGRAPHIC PROFILE



Acque del Chiampo s.p.a.
Servizio Idrico Integrato

PFOA ISOMERS: CHROMATOGRAPHIC PROFILE



Acque del Chiampo s.p.a.
Servizio Idrico Integrato

MRM TRANSITIONS: QUALIFICATION AND QUANTIFICATION

PFOA ISOMERS

413 > 369

413 > 169

+

413 > 119

413 > 219

PFOS ISOMERS

499 > 80

499 > 99

+

499 > 169


499 > 330

499 > 280



Acque del Chiampo s.p.a.
Servizio Idrico Integrato

PFOS TRANSITIONS



Acque del Chiampo s.p.a.
Servizio Idrico Integrato



PFOA TRANSITIONS

1503542-005 Smooth(Mn,4x3)
15-001504/05 89

P5-MHpA

15.95

1503542-005 Smooth(Mn,4x3)
15-001504/05 89

L-PFOA

16.10

1503542-005 Smooth(Mn,4x3)
15-001504/05 89

P6-MHpA

15.95

1503542-005 Smooth(Mn,4x3)
15-001504/05 89

P4-MHpA

16.20

15.95

F10:MRM of 4 channels,ES-
412.80 > 219.02
6.523e+005

F10:MRM of 4 channels,ES-
412.80 > 369
5.909e+006

F10:MRM of 4 channels,ES-
412.80 > 169.04
4.082e+005

F10:MRM of 4 channels,ES-
412.80 > 119.00
8.493e+004



Acque del Chiampo s.p.a.
Servizio Idrico Integrato



PFAAs



AQUEOUS SAMPLE 100%



SOLUTION?



Acque del Chiampo s.p.a.
Servizio Idrico Integrato

AMPHIPHILIC MOLECULES




Figure 1.1 Simplified surfactant structure.




Figure 1.2 Typical micelle configurations.




Acque del Chiampo s.p.a.
Servizio Idrico Integrato

PFAAs → AMPHIPHILIC MOLECULES



AGGREGATES



CONCENTRATION

PERFLUOROALKYL CHAIN LENGTH

TEMPERATURE

Na^+ CONCENTRATION (ANIONIC T.)

TURBIDITY



Acque del Chiampo s.p.a.
Servizio Idrico Integrato

SURFACE ACTIVE AGENTS




Figure 2.4 Four region isotherms of surfactant adsorption.



Acque del Chiampo s.p.a.
Servizio Idrico Integrato




Figure 2.5 Surfactant adsorption on solid surfaces.



Acque del Chiampo s.p.a.
Servizio Idrico Integrato



Use of strong anion exchange resins for the removal of perfluoroalkylated substances from contaminated drinking water in batch and continuous pilot plants

142

A. Zaggia et al. / Water Research 91 (2016) 137–146




Fig. 2. Different TEM images on microtomized sections of Purolite® AS20E saturated with PFOA. a), b) and c) show three different aggregates at the same magnification (the bar corresponds to 500 nm); d) shows aggregate c) at an additional magnification of 10× (bar corresponds to 50 nm).



Acque del Chiampo s.p.a.
Servizio Idrico Integrato

PFAAs



AMPHIPHILIC MOLECULES



AGGREGATES - ADSORPTION



**TURBIDITY – SOLIDS IN SUSPENSION
ADHESION**

**REPRESENTATIVE SAMPLE VOLUME
(ng/L!!)**



**Acque del Chiampo s.p.a.
Servizio Idrico Integrato**

PFAAs




SAMPLING?



Acque del Chiampo s.p.a.
Servizio Idrlico Integrato

CONSERVATION – PFAAs SOLUTION STABILITY



PFOS
PFNA
PFDA
PFUnA
PFDoA



con Met-OH \geq 30%
conservazione T \leq 6°C



Acque del Chiampo s.p.a.
Servizio Idrico Integrato

PFAAs

SPE



RECOVERY

ANALYTICAL
TIME

CLEAN-UP

QL

DIRECT INJECTION



SOLUTION?

MATRIX EFFECT

SOLIDS

VISCOSITY
DENSITY



Acque del Chiampo s.p.a.
Servizio Idrico Integrato

A Never-Ending Story of Per- and Polyfluoroalkyl Substances (PFASs)?

Zhanyun Wang,[†] Jamie C. DeWitt,[‡] Christopher P. Higgins,[§] and Ian T. Cousins*,^{||} 


[†]Institute for Chemical and Bioengineering, ETH Zurich, CH-8093 Zurich, Switzerland

[‡]Department of Pharmacology and Toxicology, Brody School of Medicine, East Carolina University, Greenville, North Carolina 27834, United States

[§]Department of Civil and Environmental Engineering, Colorado School of Mines, 1500 Illinois Street, Golden, Colorado 80401, United States

^{||}Department of Environmental Science and Analytical Chemistry (ACES), Stockholm University, SE-10691 Stockholm, Sweden

ABSTRACT: More than 3000 per- and polyfluoroalkyl substances (PFASs) are, or have been, on the global market, yet most research and regulation continues to focus on a limited selection of rather well-known long-chain PFASs, particularly perfluorooctanesulfonate (PFOS), perfluorooctanoic acid (PFOA) and their precursors. Continuing to overlook the vast majority of other PFASs is a major concern for society. We provide recommendations for how to proceed with research and cooperation to tackle the vast number of PFASs on the market and in the environment.



Acque del Chiampo s.p.a.
Servizio Idrico Integrato

References

World Health Organization, Keeping our water clean: the case of water contamination in Veneto region, Italy, ISBN 978 92 890 5246 7, 2017

Shoemaker J.A., Grimmet P.E., Boutin B.K., EPA 537:2009 – Determination of selected perfluorinated alkyl acids in drinking water by solid phase extraction and liquid chromatography/tandem mass spectrometry (2009)

Valsecchi S., Mazzoni M., Polesello S., Analisi multiresiduale LC-MS mediante arricchimento in linea del campione (on line SPE/UHPLC-ESI-MS/MS) per la determinazione di acidi perfluoroalchilcarbossilati e perfluoroalchilsolfonati nelle acque dolci naturali. Notiziario dei metodi analitici IRSA-CNR, vol.1 (2013)

Buck R.C. et al., Perfluoroalkyl and Polyfluoroalkyl substances in the environment: terminology, classification, and origins. Integrated environmental assessment and management, vol. 7, n. 4 (2011)

Capri S., Polesello S., Guzzella L., Carere M., Criticità connesse all'applicazione della nuova Direttiva 2013/39/UE sulle sostanze prioritarie negli ambienti acquatici. Notiziario dei metodi analitici IRSA-CNR, vol.1 (2013)



BIBLIOGRAPHY

Benskin J.P., De Silva A.O., Martin J.W., Isomer profiling of perfluorinated substances as a tool for source tracking: a review of early findings and future applications. *Reviews of environmental contamination and toxicology*, n. 208 (2010)

Esparza X. et al., Analysis of perfluorinated phosphonic acids and perfluorooctane sulfonic acid in water, sludge and sediment by LC-MS/MS. *Talanta*, n. 86 (2011)

Zaggia A., Conte L., Falletti L., Fant M., Chiorboli A., Use of strong anion exchange resins for the removal of perfluoroalkylated substances from contaminated drinking water in batch and continuous. *Water Research* 91 (2016) 137-146

Pellizzaro A., Zaggia A., Fant M., Conte L., Falletti L., Identification and quantification of linear and branched isomers of perfluorooctanoic and perfluorooctane sulfonic acids in contaminated groundwater in the Veneto region. *Journal of Chromatography A*, 1533 (2018) 143-154

Gao X., and Chorover J., Adsorption of perfluorooctanoic acid and perfluorooctanesulfonic acid to iron oxide surfaces as studied by flow-through ATR-FTIR spectroscopy. *Environ. Chem.* (2012) 9, 148–157



Acque del Chiampo s.p.a.
Servizio Idrico Integrato

BIBLIOGRAPHY

Wang Z. et al., A never-ending story of per- and polyfluoroalkyl substances (PFASs)? Environ. Sci. Technol. 2017, 51, 2508 – 2518

Wang Z. et al., Hazard assessment of fluorinated alternatives to long-chain perfluoroalkyl acids (PFAAs) and their precursors: status quo, ongoing challenges and possible solutions. Environment international 75 (2015) 172-179

Wang Z. et al., Fluorinated alternatives to long-chain perfluoroalkyl carboxylic acids (PFCAs), perfluoroalkan sulfonic acids (PFASs) and their potential precursors. Environment international 60 (2013) 242-248

Farn R. J., Chemistry and technology of surfactants (2006)

Yuan C. L. et al. Study on characteristics and harm of surfactants. J. Chem. Pharm: Res., 2014, 6(7): 2233-2237



Acque del Chiampo s.p.a.
Servizio Idrico Integrato

fantom@acquedelchiampospa.it

pellizzaroa@acquedelchiampospa.it



Acque del Chiampo s.p.a.
Servizio Idrico Integrato